

TORSO MEASUREMENT CHART



2424 E. University Drive, Phoenix, Arizona, U.S.A. 85034-6911
 800-421-5647 • 602-267-7854 • FAX: 602-273-6931
 e-mail: bio-con@bio-con.com • worldwideweb: http://www.bio-con.com

Patient _____
 Last Name First Name
 Bio-Concepts Use Only _____
 Date _____ Age _____ Sex _____

Garment color _____
 Insert _____
 Lining _____

Zipper: Front or Back Hook & loop zipper stop tab
 Regular neck
 Scoop neck: Distance below sternal notch _____
 Turtle neck: 3/4" 1" 1 1/2" 2" 2 1/2"
 Expansion panel: Left side Right side
 (default expansion panel location is lateral)
 Soft armpit panels (N/C) Waist hook & loop tabs
 Mesh armpit panels (N/C) to attach to pants (N/C)

Refer to the Bio-Concepts Measuring Manual for detailed procedures, additional instructions, and example measuring charts.

Additional instructions or comments:


CIRCLE ITEMS ORDERED

Bio-Concepts Use Only: F=

For a Vest (Items 30A, 31)
 required: A, B, C, D, G, H, b, c, d
 for adult female: brassiere cup size torso circumference over breasts
 optional: K, k (see illustrations below)


For a Body Brief or Suit (Items 24, 25, 26, 27, 28, 29)
 required: All, except K and k
 for adult female: brassiere cup size torso circumference over breasts
 optional: None

For other torso garments refer to the Measuring Manual




For Adult Females:
 Brassiere cup size: _____
 Chest circumference (over nipples): _____

CIRCUMFERENCES
 Adult Female Torso




For adult females C is directly below the breasts. Don't forget to take an additional circumference of the torso over the breasts.

CIRCUMFERENCES
 Male & Child Torso




For adult males and children, C is halfway between D and B.

HOW TO MEASURE CIRCUMFERENCES G & H




Measure H & G with arm down and place the measuring tape distal (away) from the acromion.

DISTANCES
 Adult Female Torso




b, c, d, and k are measured from Z to B, C, D, and K down the front of the torso.

DISTANCES
 Male & Child Torso



K and k are OPTIONAL. They are to extend Vests below waist line, not for Body Briefs and Body Suits.

ADDITIONAL DISTANCES & CIRCUMFERENCES
 Body Briefs & Body Suits Only



E and F are evenly spaced between the waist and the gluteal fold.

A
B
C
D
E
F
G L
H R
I L
J R
K
b
c
d
e
f
g
k