

HAND MEASUREMENT CHART

Bio-Concepts Custom Pressure Garments are available only under Physician's Order



2424 East University Drive, Phoenix, Arizona, U.S.A. 85034-6911
 VOICE: 800-421-5647 / 602-267-7854 • FAX: 800-650-9424 / 602-273-6931
 e-mail: bio-con@bio-con.com • worldwideweb: http://www.bio-con.com

PLEASE DO NOT WRITE IN MARGINS

Patient Last Name _____ First Name _____

Date _____

Garment color _____

Insert _____

Lining _____

Left hand zipper: Dorsal Palmar Ulnar Radial
 With hook & loop zipper stop tab

Right hand zipper: Dorsal Palmar Ulnar Radial
 With hook & loop zipper stop tab

Glove and Gauntlet zippers normally open at the proximal end of the garment and extend to just proximal to the knuckles (MCPs).

Finger/thumb tips closed Finger/thumb tips open


Refer to the Bio-Concepts Measuring Manual for detailed procedures, additional instructions, and example measuring charts.

Additional instructions or comments:

CIRCLE ITEMS ORDERED

Bio-Concepts Use Only:

F =



HAND TRACINGS ARE REQUIRED for all Gloves, Gauntlets, & Mittens

Use the UPPER EXTREMITY CHART

for adding sleeves to gloves & gauntlets.


Refer to the Bio-Concepts Measuring Manual for instructions on measuring contracted fingers.

REQUIRED MEASUREMENTS

Gloves: All

Mittens: All, except C thru J

Gauntlet: A1, A, B, K, d, & i




	LEFT HAND	RIGHT HAND
A1		
A		
B		
C		
D		
E		
F		
G		
H		
I		
J		
K		
a		
b		
c		
d		
e		
f		
g		
h		
i		

TO END OF FINGER TO END OF GLOVE FOR OPEN TIPS


PLEASE DO NOT WRITE IN MARGINS

REQUIRED in conjunction with HAND MEASUREMENT CHART • Position and trace patient's hand within outline

Patient _____ Date _____



If the glove is to have open finger tips, mark the desired ends of the fingers on the tracing with tic marks.



Please do not send a hand tracing without a scale

0.8x4.00

Draw a 3-inch long line in this box

PLEASE DO NOT WRITE IN MARGINS

PLEASE DO NOT WRITE IN MARGINS

Take measurements a-i on the Hand Measurement Chart from the hand tracing.


Place wrist crease on this line

REQUIRED in conjunction with HAND MEASUREMENT CHART • Position and trace patient's hand within outline

Patient _____ Date _____

0.8x4.00


Correct



Incorrect




If the glove is to have open finger tips, mark the desired ends of the fingers on the tracing with tic marks.



Draw a 3-inch long line in this box

PLEASE DO NOT WRITE IN MARGINS
Please do not send a hand tracing without a scale



Take measurements
a-i on the Hand
Measurement Chart
from the hand tracing.

Place wrist crease on this line

PLEASE DO NOT WRITE IN MARGINS